



HOW TO SAVE A DEBUG LOG _____

How to Save a Debug Log with RADAR I and RADAR II:

- 1 Use the **MENU/PREV**, the **↑** and **↓** keys and the **ENTER** key to go to the **MAIN MENU / DIAGNOSTICS MENU / SAVE DEBUG** dialogue and press the **ENTER** key.
- 2 Insert a floppy into the disk drive. At the **SAVE DEBUG INFO** prompt type in a name for the debug file or leave it as the default **DEBUG.OUT** and press the **ENTER** key. Insert the floppy into a Mac/PC, and email the relevant files to support@izcorp.com with a detailed account of the issue.